



Ready to wear  
compression  
garments.

[readyfitgarments.com.au](https://readyfitgarments.com.au)

**readyfit**  
By Second Skin

# Readyfit compression garments



## Compression gloves

Designed to offer mild compression to the hand, the Readyfit compression glove provides both durability and comfort.

Optional zipper provides easy application and reduced friction on the skin.



### Features:

- Soft gussets to reduce irritation in the finger and web spaces
- Offered with or without zipper
- Can be worn in conjunction with a Readyfit arm sleeve
- Offered in light, dark and black colour options and a range of sizes

### Effective for:

- Managing mild oedema
- Improving the appearance of scarring
- Managing stable or healed skin grafts
- Managing partial thickness burns



TGA Approved

## Compression arm sleeves

Extends from the axilla to the wrist providing effective skin contact up the arm.

Optional zipper provides easy application and reduced friction on the skin.



### Features:

- Made from light open weave fabric for increased cooling properties
- Includes soft fabric in the front of the elbow for increased comfort
- Offered with or without zipper
- Can be worn in conjunction with a Readyfit glove
- Offered in light, dark and black colour options and a range of sizes

### Effective for:

- Managing mild oedema through the arm
- Improving the appearance of scarring
- Managing stable or healed skin grafts
- Managing partial thickness burns



TGA Approved



Order online: [readyfitgarments.com.au](https://www.readyfitgarments.com.au)



## Knee high compression socks

Designed with vertical seams to improve the fit and contour over the ankle to maximise skin contact over scarring and facilitate the flow of fluid out of the foot, ankle and calf.

Optional zipper provides easy application and reduced friction on the skin.



### Features:

- Made from light open weave fabric for increased cooling properties and comfort
- Offered in open or closed toe options
- Offered with or without zipper
- Offered in light, dark and black colour options and a range of sizes

### Effective for:

- Managing mild lower limb oedema
- Improving the appearance of scarring
- Fragile or sensitive skin



TGA Approved

## Ankle high compression socks

Shaped to ensure effective skin contact and compression and a great option when a full-length sock is not required.

### Features:

- Made from light open weave fabric for increased cooling properties and comfort
- Offered in open or closed toe options
- Offered in light, dark and black colour options and a range of sizes

### Effective for:

- Localised oedema in the foot
- Managing scarring over the dorsum of the foot
- Fragile or sensitive skin



TGA Approved

# Readyfit reusable face masks



## Soft fabric face mask

Washable and reusable face mask, shaped to contour the face. Developed from the same soft fabric used within our custom-made garments to provide a comfortable feel and durability.

### Features:

- Durable and machine washable for up to 3 months
- Double polyester and elastane layer
- Contoured shaping to ensure good contact around the face
- Soft fabric finish does not imprint on the skin
- Comfortable ear loops to ensure the mask stays in place
- Available in adult and small adult/child sizes in light and black colours

*Disclaimer: This is not a surgical grade/N95 respirator mask. It is intended for general use. Not to be worn for more than 4 hours at a time.*



## Breathable face mask

Washable and reusable face mask made from breathable polypropylene layers, which allow for longer wear duration.

### Features:

- Double layer of 100% polypropylene fabric
- Maintains effectiveness after washing (approximately 10 washes)
- Contoured shaping to ensure good contact around the face
- Comfortable ear loops to ensure the mask stays in place
- Water resistant until washed
- Available in adult and small adult/child sizes in a mix of colours

*Disclaimer: This is not a surgical grade/N95 respirator mask. It is intended for general use. Level 3 rating for differential pressure testing AS 4381:2015 maintained after washing.*



## Did you know?

Our fabrics are treated with a patented dynamic cooling technology, **HeiQ Smart Temp**. Designed to keep clients cooler, dryer and more comfortable.

For over 30 years we have been a leader in the design and manufacture of compression garments.



We have an established reputation for quality and we work closely with health professionals and clients to ensure we provide the solution that will work for them.



Our experienced team consists of occupational therapists, physiotherapists, clinic management, design, production and administrative staff who are driven to improve the daily lives for all our clients. Together we provide garments that are comfortable, durable, cost effective and complemented by valuable expertise and assistance.



## Why choose Readyfit?



**Consistent effective medical compression**



**Accurate fitting**



**Leading edge design and manufacture**



**Features and options**



**High quality fabrics**



**Support and assistance**

Our garments are made with the needs of our clients in mind.



If you require a custom fit solution, or fall between a size in the Readyfit range, we offer customised garments to deliver your required outcome.

We have been designing and manufacturing customised garments for over 30 years and our skilled team would be happy to talk about your specific clinical needs and goals.

Visit: [secondskin.com.au](https://secondskin.com.au) for more information.



**T +61 8 9201 9455**

**E [readyfit@secondskin.com.au](mailto:readyfit@secondskin.com.au)**

**[readyfitgarments.com.au](https://readyfitgarments.com.au)**

Unit 15, 386 Scarborough Beach Road  
Osborne Park WA 6017

Follow us 



**Made in Australia**

Second Skin is a registered National Disability Insurance Scheme (NDIS) provider.



**Second Skin**  
CUSTOMISED COMPRESSION GARMENTS & SPLINTS



Second Skin is a leader in the design, manufacture and service delivery of custom made medical compression garments and dynamic splints.