



Medisafe

AQ Box Cloud Temperature Monitoring



Why choose us

Designed and proudly built in Australia since 2010, the Medisafe Plus family of medical refrigeration units and freezers are purpose-built for vaccine storage and follow the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines.

All Medisafe Plus vaccine storage units come complete with the built-in Medisafe Plus cloud monitoring system and 12 months' free access to the portal. Known as the AQ Box, the system is an online temperature monitoring and alerting device which automatically records and sends temperature data to the cloud for easy viewing and reporting via the AQ Box portal. Users log in with their unique username and password to easily view and download or print reports for compliance requirements.

The AQ Box comes standard with all Medisafe Plus refrigeration units or can be fitted to your existing vaccine storage unit.

Dual Temperature Probes

The AQ Box is unlike other monitoring devices on the market which measure the cabinet temperature only. The AQ Box is uniquely equipped with two temperature probes providing accurate measurements of both the cabinet temperature and simulated product temperature.

Remote Access

Temperature data can be accessed via the AQ Box portal on any web-enabled smartphone, tablet or computer. Each unit comes with 12 months' access to the cloud portal*. Temperature breaches and power failure notifications are sent in real time via text and email.

*Averm Quirks retains ownership of the AQ Box and will manage the service agreement with the telecommunications provider.

* Charges are applicable after the initial 12 month term.

The AQ Box is the self-managed temperature monitoring solution that automatically logs data to the cloud.


Downloadable reports are available with nominated date ranges for ease of compliance with regulatory requirements and can include the duration of any temperature breaches.

Power is supplied by a 10amp power source. However, if power is disrupted, the built-in backup battery will ensure uninterrupted temperature recording of up to 7 hours* during a power outage.

*Back up times are based on a new unit and will reduce with usage.

The AQ Box ensures that expensive vaccines are protected, your data is secure and that you are meeting the national Australian Department of Health's 'Strive for 5' guidelines.

Temperature Monitoring Device Working




AVEM Field Test Unit 2:
9 Business Park Drive, Lynbrook, Victoria, Australia

Device Working:
Device in now working perfectly

Report Time:
2021-08-24 12:58:40 +1000 AEST

View Details

Email notification of device status with time and date



Product and fridge temperature graph

Features

- Dual RTD Class A temperature probes for cabinet and simulated vaccine temperatures
- Accuracy $\pm 0.15^{\circ}\text{C}$ @ 0°C , $\pm 0.35^{\circ}\text{C}$ @ 100°C
- 7-hour rechargeable battery back up
- Cache of data if there is a connectivity issue
- Full web reporting and easily accessible data
- Minimum and Maximum reporting available
- Data accessible up to 90 days

Optional Extras

- Wi-Fi extender (for connection to Wi-Fi network)
- Ethernet connection

DIMENSIONS (HxWxD)	110 x 110 x 32mm
POWER CONSUMPTION /24 HOURS	1.2 kWh
CONNECTIVITY	4G

Local distributor:

